

FGF2 DISC™ - Standard

Product Information Sheet




StemCultures | ☎ +1 518 621 0848 | ✉ support@stemcultures.com | 📍 1 Discovery Drive | Rensselaer, New York 12144 | United States
🌐 www.stemcultures.com | 🐦 @StemCultures | 📺 Stem Cultures

Product Description

FGF2 DISCs™ are inert, non-degradable, biocompatible hydrogels that release defined levels of native FGF2 stably into culture medium over the course of one week. FGF2 DISCs™ are loaded with StemBeads® FGF2, PLGA microbeads that encapsulate native, human FGF2 as controlled-release technology. FGF2 DISCs™ are easy to add and remove, giving scientists enhanced control of growth factor levels in their cultures. Controlled delivery and stable levels overcome the 4.5 hour half-life (Figure 1) of FGF2 and improve cell cultures while saving researchers valuable time and resources.

FGF2 DISCs™ have been tested in medium such as mTeSR™, mTeSR™ Plus, StemFlex™, Essential 8™, and neuronal medium with enhanced cellular profiles. FGF2 DISCs™ can be combined with other StemCultures products.

Figure 1




Product Information

Catalog #	Product Name	Storage	Expiration	DISC™ Size	Recommended Well/Plate Size
DSC500S, DSC500-48	FGF2 DISC™ (Standard)	4°C	6 months, as specified on label	2-3 mm diameter, dry 5-6 mm diameter, rehydrated	6 well, 12 well

Suggested Protocol for Use in Pluripotent Stem Cell (PSC) Cultures

- Using aseptic cell culture technique, wipe DISC™ container with 70% ethanol and place into a biosafety cabinet before opening.
- Passage cells and add culture medium to wells.
- Using sterile forceps, transfer each DISC™ into a culture well containing 2 mL of medium (see Release Data section on page 2 for more information on release and adjustments).
Note: As DISCs™ rehydrate, they will swell and become transparent. Embedded StemBeads® will be visible under a microscope. An image of a DISC™ under a microscope is on page 2 for reference.
- Every 2-3 days, replace only the medium (use a low powered vacuum or a pipette), leaving the original DISC™ in the well.
- After about 7 days total, passage cells into a new culture dish and add a new DISC™. The old DISC™ can be removed using a low powered vacuum + pipette tip.

Recommended PSC Culture Schedule



Note: Different PSC lines, culture densities, and media may require adjusted schedules.

Please reach out to support@stemcultures.com for ordering and technical support.

Release Data

StemCultures offers two product sizes to accommodate standard culture plate sizes and to allow more flexibility with the DISC™ release. We recommend a release of 10 ng/mL when culturing iPSCs. However, to fit other needs, the release can be adjusted slightly based on the amount of medium and the number of DISCs™ that are added. See the chart below for reference.

DISC™ Size	Volume of Medium Added	Number of DISCs™ per well	Release in Volume of Medium Added	Example Plate Size
Standard (DSC500)	2 mL	1	10 ng/mL	6 well
Standard (DSC500)	1 mL	1	20 ng/mL	12 or 24 well
Standard (DSC500)	0.5 mL	1	40 ng/mL	24 well
Standard (DSC500)	2 mL	2	20 ng/mL	6 well
Standard (DSC500)	2 mL	3	30 ng/mL	6 well
Standard (DSC500)	2 mL	4	40 ng/mL	6 well

Visual of DISCs™ in Culture

FGF2 DISCs™ will become transparent when rehydrated. Cells under the DISC™ and the embedded StemBeads® will be visible under a microscope as seen in Figure 2 below.

Figure 2

